

pTriEx-4 Neo Vector

Baculovirus Locus	polh
Promoters	CMV immediate early p10 T7/lac
N-terminal fusion options	His•Tag S•Tag
C-terminal fusion options	HSV•Tag His•Tag
Cloning options	polylinker
Selectable marker (mammalian cells)	Neomycin ^R




pTriEx-4 Neo sequence landmarks

CMV ie enhancer/promoter	1021–1597
Vertebrate transcription start	1598
T7 promoter	1931–1947
T7 transcription start	1948
lac operator	1952–1972
p10 promoter region	1986–2099
p10 transcription start	2030–2031
Multiple cloning sites (Nco I–Bsu36 I)	2112–2456
His•Tag [®] coding sequence	2120–2137
S•Tag [™] coding sequence	2147–2191
HSV•Tag [®] coding sequence	2378–2413
His•Tag [®] coding sequence	2420–2443
CITE sequence	2505–3004
Rabbit globin terminator region	3861–4103
T7 terminator	4104–4151
pUC origin	5080
bla coding sequence	5679–6539

The pTriEx[™]-4 Neo vector¹ (Cat. No. 70933-3) is designed to allow rapid characterization of target genes in *E. coli*, insect and vertebrate cells, and to allow rapid selection of stable transfected vertebrate cells expressing high levels of target gene. Expression in vertebrate cells is mediated by the CMV immediate early enhancer and promoter². The drug selection marker is expressed under the control of the EMC virus Cap-Independent Translation Enhancer (CITE) sequence (or IRES), allowing rapid selection of transfected vertebrate cells using the drug G 418 or neomycin sulfate. For expression in insect cells, pTriEx-4 Neo contains flanking baculovirus sequences to permit the generation of recombinant baculoviruses using the BacVector[®] System. In baculovirus-infected insect cells, expression is driven by the very late p10 promoter. Expression in *E. coli* is regulated by the tightly controlled T7/lac promoter. Expression can be induced in hosts such as NovaBlue by infecting with λCE6, a phage that constitutively expresses T7 RNA polymerase from the λ_{p_L} and λ_{p_T} promoters. Alternatively, pTriEx recombinant plasmids can be transferred into a (DE3)pLacI host that allows induction with IPTG. Native protein can be expressed by cloning into the *Nco* I site, or generated by cloning into the *Psh* A I or *Sma* I sites and cleaving the fusion protein with enterokinase or thrombin, respectively.

¹ patent pending

² The CMV promoter is covered by U.S. Patent nos. 5,168,062 and 5,385,839 issued to the University of Iowa Research Foundation and is licensed for research use only.



pTriEx-4 Neo Restriction Sites

Enzyme	# Sites	Locations	Enzyme	# Sites	Locations	Enzyme	# Sites	Locations		
AatII	6	1140 1193 1276 1462 1583 4609	DpnI	29		PvuII	1	2362		
AccI	3	244 2328 2366	DraI	4	429 4048 4418 6329	RcaI	1	3928		
AcI	61		DraIII	2	2446 2794	RsaI	18			
AflIII	7	399 2369 2746 2921 4679 4829 5018	DrdI	2	3204 5125	RsrII	2	1921 3693		
AhdI	2	499 5752	DsaI	3	1743 2112 2965	SacI	1	2289		
AluI	21		EaeI	9	1011 2352 2734 3047 3083 3257 3648 3675 6140	SacII	1	1746		
Alw26I	8	285 949 1449 1659 1897 4014 5813 6589	EagI	2	2352 3083	Sall	2	243 2327		
AlwI	16		EarI	6	51 547 2143 3521 3731 6547	SapI	2	3521 3731		
AlwNI	3	1865 2231 5433	Ecl136II	1	2287	Sau3AI	29			
ApaI	2	1712 2548	Eco47III	1	4286	Sau96I	19			
ApaLI	3	2908 5331 6419	Eco57I	5	2876 3323 3755 5565 6419	Scal	1	6232		
ApoI	12		EcoO109I	6	1709 1989 2544 2900 2954 4120	ScrFI	23			
AscI	1	2315	EcoRI	1	2299	SfaNI	16			
AvaI	4	1732 2214 2405 2414	EcoRII	11		Sfcl	9	1857 1943 2321 2679 3029 3226 5282 5473 5993		
Avall	7	1921 1989 2209 3693 4237 5890 6112	EcoRV	1	2280	Sgfl	1	4874		
AvrII	1	2582	EheI	1	3178	SmaI	1	2216		
BamHI	1	2293	FauI	9	1108 1134 1301 1529 1712 1738 1775 3024 3300	SnaBI	2	1355 4328		
BanI	7	1480 2333 2724 2872 3176 3211 5700	Fnu4HI	45		SphI	1	4162		
BanII	4	1712 2289 2548 3542	FokI	6	2873 3501 3526 5718 5899 6186	Sse8387I	1	2325		
BbsI	5	498 2406 2557 2671 4787	FseI	1	1017	Sspl	4	425 4422 4619 6556		
BbvI	25		FspI	3	659 3279 5974	Styl	3	2112 2582 4115		
Bcgl	1	6291	HaeII	3	3180 4288 5265	Swal	1	4418		
Bfal	12		HaeIII	25		Tail	25			
BglI	5	1105 1227 1298 2733 5872	Hgal	9	146 503 966 1541 4228 4388 4814 5128 6278	TaqI	18			
BglII	1	2307	Hhal	24		TfiI	5	446 3662 3796 4291 4993		
Bpml	1	5822	HincII	3	245 1531 2329	Thal	16			
Bpu10I	1	416	HindIII	3	2345 2655 3003	Tsel	25			
BsaAI	5	1355 2372 2747 3481 4328	HinfI	15		Tsp45I	5	1802 3297 3603 6008 6219		
BsaHI	10	495 1137 1190 1273 1459 1580 3177 4606 4806 6289	HphI	11	183 879 2117 2420 2426 2821 3355 5822 6238 6444	Tsp509I	42			
Bsal	2	285 5813	KpnI	2	2337 2876	TspRI	14			
BsaJI	17		MaeIII	15		Tth111I	1	3295		
BsaWI	7	833 2336 2502 3208 5223 5370 6043	MbolI	21		VspI	5	1022 1930 2057 4813 5924		
BseRI	1	2260	MluI	1	399	XbaI	1	1815		
BsgI	2	1689 2482	MnII	35		XcmI	1	2138		
BsiEI	7	2355 3086 4874 4934 5357 6122 6271	MscI	1	3259	XhoI	1	2414		
BsiHKAI	8	2289 2421 2912 3290 3480 5335 6338 6423	MseI	38		XmnI	2	2643 6351		
BsII	19		MslI	8	950 1380 1837 3614 4057 6004 6163 6522	Enzymes that do not cut pTriEx-4 Neo:				
BsmBI	1	949	MspA1I	8	655 1745 2156 2362 2701 5359 5604 6387	AflIII	BclI	Bpu1102I	BsaBI	BspEI
BsmFI	11		MspI	31		BstEII	EcoNI	HpaI	NheI	NruI
BsmI	2	2577 2616	MunI	2	4317 4662	PmeI	SanDI	SexAI	SfiI	SgrAI
Bsp1286I	13		MwoI	29		SpeI	SrfI	StuI	SunI	
BspLU11I	2	2921 5018	NarI	1	3177					
BspMI	6	1686 2216 2314 2770 3064 3445	NciI	12						
BsrBI	4	1728 1960 3790 4951	NcoI	1	2112					
BsrDI	6	79 2542 3410 3979 5813 5987	NdeI	4	1249 4028 4088 4096					
BsrFI	8	42 781 1013 1806 2336 3496 3677 5832	NgoAIV	3	781 1013 3677					
BsrGI	3	49 768 4658	NlaIII	20						
BsrI	16		NlaIV	21						
BssHII	1	2315	NotI	1	2352					
BssSI	4	2288 3769 5190 6416	Nsil	1	184					
Bst1107I	1	2367	NspI	4	2916 2925 4162 5022					
BstXI	2	167 3868	NspV	1	2342					
BstYI	9	2293 2307 2410 2888 3348 5658 5669 6377 6394	Pacl	1	2061					
Bsu36I	1	2456	PfiMI	2	2182 2884					
Cac8I	29		PinAI	1	2336					
Clal	1	3892	PleI	10	150 1419 1931 2180 2265 2825 3600 4661 5396 5741					
CviJI	101		PmlI	2	2372 2747					
Ddel	10	416 2283 2304 2443 2456 3840 3886 5292 5709 6249	PshAI	1	2256					
			Psp1406I	2	5978 6351					
			Psp5II	1	1989					
			PstI	3	2325 3033 3230					
			PvuI	2	4874 6122					