

# pVAX1<sup>©</sup>

Catalog no. V260-20

**Version B**

*010124  
25-0256*



**U.S. Headquarters:**

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008  
Toll Free Tel: (800) 955-6288  
Tel: (760) 603-7200  
Fax: (760) 602-6500  
E-mail: [tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)  
Web Site: [www.invitrogen.com](http://www.invitrogen.com)

**European Headquarters:**

Invitrogen BV  
PO Box 2312, 9704 CH Groningen  
The Netherlands  
Toll Free Tel: 00800 5345 5345  
Toll Free Fax: 00800 7890 7890  
Tel: +31 (0) 50 5299 299  
Fax: +31 (0) 50 5299 281  
E-mail: [tech\\_service@invitrogen.nl](mailto:tech_service@invitrogen.nl)



# Table of Contents

TABLE OF CONTENTS .....	III
GENERAL INFORMATION .....	V
PURCHASER NOTIFICATION.....	VI
METHODS .....	1
OVERVIEW .....	1
CLONING INTO PVAX1 <sup>®</sup> .....	2
TRANSIENT TRANSFECTION .....	4
APPENDIX .....	5
PVAX1 <sup>®</sup> VECTOR.....	5
PVAX1/LACZ <sup>®</sup> MAP .....	7
TECHNICAL SERVICE .....	8
REFERENCES.....	10



## General Information

---

**Contents** 20 µg each of pVAX1<sup>®</sup> and pVAX1/*lacZ*<sup>®</sup>, lyophilized in TE, pH 8.0

---

**Shipping/Storage** Lyophilized plasmids are shipped at room temperature and should be stored at -20°C.

---

**Product Qualification** The pVAX1 and pVAX1/*lacZ* vectors are qualified by restriction digest. Restriction digests must demonstrate the correct banding pattern when electrophoresed on an agarose gel. The table below lists the restriction enzymes used to digest each vector and the expected fragments.

Vector	Restriction Enzyme	Expected Fragments (bp)
pVAX1	<i>EcoR</i> V	2999
	<i>Hind</i> III	2999
	<i>Nco</i> I	1375, 1624
	<i>Pml</i> I	2999
pVAX1/ <i>lacZ</i>	<i>EcoR</i> V	6050
	<i>Hind</i> III	6050
	<i>Nco</i> I	1624, 4426
	<i>Pml</i> I	6050

---

### Technical Service



For Technical Service, please call, write, fax or E-mail:

**U.S. Headquarters:**

Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008  
Toll Free Tel: (800) 955-6288  
Tel: (760) 603-7200  
Fax: (760) 602-6500  
E-mail: [tech\\_service@invitrogen.com](mailto:tech_service@invitrogen.com)  
Web Site: [www.invitrogen.com](http://www.invitrogen.com)

**European Headquarters:**

Invitrogen BV  
PO Box 2312, 9704 CH Groningen  
The Netherlands  
Toll Free Tel: 00800 5345 5345  
Toll Free Fax: 00800 7890 7890  
Tel: +31 (0) 50 5299 299  
Fax: +31 (0) 50 5299 281  
E-mail: [tech\\_service@invitrogen.nl](mailto:tech_service@invitrogen.nl)

---

# Purchaser Notification

---

## Introduction

Use of the pVAX1<sup>®</sup> vector is covered under a number of different licenses including those detailed below.

---

## CMV Promoter License

Use of the CMV promoter is covered under U.S. Patent Nos. 5,168,062 and 5,385,839 owned and licensed by the University of Iowa Research Foundation and may be used for research purposes only. Commercial users must obtain a license to these patents directly from the University of Iowa Research Foundation. Inquiries for commercial use should be directed to:

Brenda Akins  
University of Iowa Research Foundation (UIRF)  
214 Technology Innovation Center  
Iowa City , IA 52242  
Tel: 319-335-4549

---

## BGH PolyA Sequence License

The bovine growth hormone (BGH) polyadenylation sequence is licensed under U.S. Patent No. 5,122,458 for research purposes only. "Research purposes" means uses directed to the identification of useful recombinant proteins and the investigation of the recombinant expression of proteins, which uses shall in no event include any of the following:

- a. any use in humans of a CLAIMED DNA or CLAIMED CELL;
- b. any use in human of protein or other substance expressed or made at any stage of its production with the use of a CLAIMED DNA or a CLAIMED CELL;
- c. any use in which a CLAIMED DNA or CLAIMED CELL would be sold or transferred to another party other than Invitrogen, its AFFILIATE, or its SUBLICENSEE;
- d. any use in connection with the expression or production of a product intended for sale or commercial use; or
- e. any use for drug screening or drug development.

Inquiries for commercial use should be directed to:

Bennett Cohen, Ph.D.  
Research Corporation Technologies  
101 North Wilmot Road, Suite 600  
Tucson, AZ 85711-3335  
Tel: 1-520-748-4443  
Fax: 1-520-748-0025

---

# Methods

## Overview

### Introduction

---

pVAX1<sup>®</sup> is a 3.0 kb plasmid vector designed for use in the development of DNA vaccines. The vector was constructed to be consistent with the Food and Drug Administration (FDA) document, “Points to Consider on Plasmid DNA Vaccines for Preventive Infectious Disease Indications”, published December 22, 1996 (see **FDA “Points to Consider”** below). Features of the vector allow high-copy number replication in *E. coli* and high-level transient expression of the protein of interest in most mammalian cells (see pages 5-6 for more information). The vector contains the following elements:

- Human cytomegalovirus immediate-early (CMV) promoter for high-level expression in a wide range of mammalian cells
- Bovine growth hormone (BGH) polyadenylation signal for efficient transcription termination and polyadenylation of mRNA
- Kanamycin resistance gene for selection in *E. coli*

The control plasmid, pVAX1/*lacZ*<sup>®</sup>, is included for use as a positive control for transfection and expression in the cell line of choice.

---

### Experimental Outline

Use the following outline to clone and express your gene of interest in pVAX1.

- Consult the multiple cloning site described on page 3 to design a strategy to clone your gene into pVAX1.
  - Ligate your insert into the appropriate vector and transform into *E. coli*. Select transformants on LB plates containing 50 µg/ml kanamycin.
  - Analyze your transformants for the presence of insert by restriction digestion.
  - Select a transformant with the correct restriction pattern and use sequencing to confirm that your gene is cloned in the proper orientation.
  - Transfect your construct into the mammalian cell line of choice and test for transient expression of your protein of interest.
- 

### FDA “Points to Consider”

pVAX1 was constructed by modifying the vector, pcDNA3.1<sup>®</sup>, to accommodate the following considerations put forth by the FDA Center for Biologics Evaluation and Research (CBER) in its document, “Points to Consider on Plasmid DNA Vaccines for Preventive Infectious Diseases Indications” (Docket no. 96N-0400).

- Sequences not necessary for replication in *E. coli* or for expression of recombinant protein in mammalian cells were removed to limit DNA sequences with possible homology to the human genome and to minimize the possibility of chromosomal integration.
- The kanamycin resistance gene was substituted for the ampicillin resistance gene because aminoglycoside antibiotics are less likely to elicit allergic responses in humans.

For further details, you may download the “Points to Consider” document from the World Wide Web at the following address:

<http://www.fda.gov/cber/points.htm>

---

# Cloning into pVAX1<sup>®</sup>

---

## Introduction

A diagram is provided on the next page to help you design a cloning strategy for ligating your gene of interest into pVAX1. General considerations for transformation are listed below.

---

## General Molecular Biology Techniques

For help with DNA ligation, *E. coli* transformation, restriction enzyme analysis, DNA sequencing, and DNA biochemistry, please refer to *Molecular Cloning: A Laboratory Manual* (Sambrook *et al.*, 1989) or *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

---

## *E. coli* Strain for Transformation

Many *E. coli* strains are suitable for the propagation of this vector, including TOP10 (Catalog no. C610-00) or DH5a<sup>™</sup>-T1<sup>R</sup>. We recommend that you propagate vectors containing inserts in *E. coli* strains that are recombination deficient (*recA*<sup>1</sup>) and endonuclease A deficient (*endA*).

For your convenience, TOP10 and DH5a<sup>™</sup>-T1<sup>R</sup> *E. coli* are available as chemically competent or electrocompetent (TOP10 only) cells in a One Shot<sup>®</sup> format from Invitrogen.

Item	Quantity	Catalog no.
One Shot <sup>®</sup> TOP10 (chemically competent cells)	21 x 50 µl	C4040-03
One Shot <sup>®</sup> TOP10 Electrocomp (electrocompetent cells)	21 x 50 µl	C4040-52
One Shot <sup>®</sup> DH5a <sup>™</sup> -T1 <sup>R</sup> MAX Efficiency <sup>®</sup> (chemically competent cells)	21 x 50 µl	12297-016

---

## Transformation Method

You may use any method of your choice for transformation. Chemical transformation is the most convenient for most researchers. Electroporation is the most efficient and the method of choice for large plasmids.

---

## Maintenance of pVAX1

To propagate and maintain the pVAX1 plasmid, resuspend the vector in 20 µl sterile water to prepare a 1 µg/µl stock solution. Store the stock solution at -20°C.

Use this stock solution to transform a *recA*<sup>1</sup>, *endA* *E. coli* strain like TOP10, DH5a<sup>™</sup>-T1<sup>R</sup>, or equivalent. Select transformants on LB plates containing 50 µg/ml kanamycin. Be sure to prepare a glycerol stock of your plasmid-containing *E. coli* strain for long-term storage (see page 3).

---

## Cloning Considerations

pVAX1 is a nonfusion vector. Your insert must contain a Kozak translation initiation sequence and an initiation codon (ATG) for proper initiation of translation (Kozak, 1987; Kozak, 1991; Kozak, 1990). An example of a Kozak consensus sequence is provided below. Please note that other sequences are possible (see references above), but the G or A at position -3 and the G at position +4 are the most critical (shown in bold). The ATG initiation codon is shown underlined.

(G/A)NNATGG

Your insert must also contain a stop codon for proper termination of your gene. Please note that the *Xba* I site contains an internal stop codon (TCTAGA).

---


*continued on next page*



# Cloning into pVAX1<sup>®</sup>, continued

## Multiple Cloning Site of pVAX1

Below is the multiple cloning site for pVAX1. Restriction sites are labeled to indicate the cleavage site. The multiple cloning site has been confirmed by sequencing and functional testing. **The complete sequence may be downloaded from our web site (<http://www.invitrogen.com>)** or from Technical Service (see page 8).



\*Please note that there are two *Pme* I sites and two *BstX* I sites in the polylinker.

## Transformation of Ligation Mixture

Transform your ligation mixture into a competent *recA*<sup>1</sup>, *endA* *E. coli* strain (e.g. TOP10, DH5 $\alpha$ ) and select on LB plates containing 50  $\mu$ g/ml kanamycin. Select 10-20 clones and analyze for the presence and orientation of your insert.



We recommend that you sequence your construct with the T7 and BGH Reverse primers (Catalog nos. N560-02 and N575-02, respectively) to confirm that your gene is cloned in the proper orientation for expression and that it contains an ATG and a stop codon. See the diagram above for the sequences and location of the priming sites. The primers are available separately from Invitrogen in 2  $\mu$ g aliquots.

## Preparing a Glycerol Stock

Once you have identified the correct clone, purify the colony and make a glycerol stock for long-term storage. It is also a good idea to keep a DNA stock of your plasmid at -20°C.

- Streak the original colony out for single colonies on an LB plate containing 50  $\mu$ g/ml kanamycin. Incubate the plate at 37°C overnight.
- Isolate a single colony and inoculate into 1-2 ml of LB containing 50  $\mu$ g/ml kanamycin.
- Grow the culture to mid-log phase (OD<sub>600</sub> = 0.5-0.7).
- Mix 0.85 ml of culture with 0.15 ml of sterile glycerol and transfer to a cryovial.
- Store at -80°C.

# Transient Transfection

---

## Introduction

Once you have verified that your gene is cloned in the correct orientation and contains an initiation ATG and a stop codon, then you are ready to transiently transfect your mammalian cell line of choice to check for protein expression. We recommend that you include the positive control vector and a mock transfection (negative control) to evaluate your results.

---

## Plasmid Preparation

Plasmid DNA for transfection into eukaryotic cells must be very clean and free from phenol and sodium chloride. Contaminants will kill the cells, and salt will interfere with lipids, decreasing transfection efficiency. We recommend isolating plasmid DNA using the S.N.A.P.<sup>™</sup> MiniPrep Kit (10-15  $\mu$ g DNA, Catalog no. K1900-01), the S.N.A.P.<sup>™</sup> MidiPrep Kit (10-200  $\mu$ g DNA, Catalog no. K1910-01), or CsCl gradient centrifugation.

---

## Methods of Transfection

For established cell lines (e.g. HeLa), please consult original references or the supplier of your cell line for the optimal method of transfection. Follow exactly the protocol for your cell line. Pay particular attention to medium requirements, when to pass the cells, and at what dilution to split the cells. Further information is provided in *Current Protocols in Molecular Biology* (Ausubel *et al.*, 1994).

Methods for transfection include calcium phosphate (Chen and Okayama, 1987; Wigler *et al.*, 1977), lipid-mediated (Felgner *et al.*, 1987; Felgner *et al.*, 1989), and electroporation (Chu *et al.*, 1987; Shigekawa and Dower, 1988). Invitrogen offers the Calcium Phosphate Transfection Kit (Catalog no. K2780-01) and a large selection of reagents for transfection. For more information on the reagents available, please visit our World Wide Web site ([www.invitrogen.com](http://www.invitrogen.com)) or call Technical Service (see page 8).

---

## Positive Control

pVAX1/*lacZ* is provided as a positive control vector for mammalian transfection and expression (see page 7). It may be used to optimize transfection conditions for your cell line. The gene encoding  $\beta$ -galactosidase is expressed in mammalian cells under the control of the CMV promoter. A successful transfection will result in  $\beta$ -galactosidase expression that can be easily assayed (see below).

---

## Assay for $\beta$ -galactosidase Activity

You may assay for  $\beta$ -galactosidase expression by activity assay using cell-free lysates (Miller, 1972) or by staining the cells for activity. Invitrogen offers the  $\beta$ -Gal Assay Kit (Catalog no. K1455-01) and the  $\beta$ -Gal Staining Kit (Catalog no. K1465-01) for fast and easy detection of  $\beta$ -galactosidase expression.


---

## Appendix

### pVAX1<sup>®</sup> Vector

#### Map of pVAX1

The figure below summarizes the features of the pVAX1 vector. The sequence for pVAX1 is available for downloading from our World Wide Web site (<http://www.invitrogen.com>) or from Technical Service (see page 8).



#### Comments for pVAX1<sup>®</sup>: 2999 bp

CMV promoter: bases 137-724  
T7 promoter/priming site: bases 664-683  
Multiple cloning site: bases 696-811  
BGH reverse priming site: bases 823-840  
BGH polyadenylation signal: bases 829-1053  
Kanamycin resistance gene: bases 1226-2020  
pUC origin: bases 2320-2993

*continued on next page*

## pVAX1<sup>®</sup> Vector, continued

### Features of pVAX1

pVAX1 (2999 bp) contains the following elements. All features have been functionally tested.

<b>Feature</b>	<b>Benefit</b>
Human cytomegalovirus (CMV) immediate-early promoter/enhancer	Permits efficient, high-level expression of your recombinant protein (Andersson <i>et al.</i> , 1989; Boshart <i>et al.</i> , 1985; Nelson <i>et al.</i> , 1987)
T7 promoter/priming site	Allows for <i>in vitro</i> transcription in the sense orientation and sequencing through the insert
Multiple cloning site	Allows insertion of your gene and facilitates cloning
BGH reverse priming site	Permits sequencing through the insert
Bovine growth hormone (BGH) polyadenylation signal	Efficient transcription termination and polyadenylation of mRNA (Goodwin and Rottman, 1992)
Kanamycin resistance gene	Selection of vector in <i>E. coli</i> (Davies and Smith, 1978)
pUC origin	High-copy number replication and growth in <i>E. coli</i> (Bolivar <i>et al.</i> , 1977; Bolivar <i>et al.</i> , 1977)


# pVAX1/*lacZ*<sup>®</sup> Map

## Description

pVAX1/*lacZ* is a 6050 bp control vector containing the gene for  $\beta$ -galactosidase. The vector was constructed by cloning a 3.1 kb *Pst* I-*Xba* I fragment containing the *lacZ* gene into the *Pst* I-*Xba* I site of pVAX1.

## Map of Control Vector

The figure below summarizes the features of the pVAX1/*lacZ* vector. **The complete nucleotide sequence for pVAX1/*lacZ* is available for downloading from our World Wide Web site (<http://www.invitrogen.com>) or by contacting Technical Service.** See the next page for more information.



### Comments for pVAX1/*lacZ*<sup>®</sup>: 6050 bp

CMV promoter: bases 137-724

T7 promoter/priming site: bases 664-683

LacZ ORF: bases 773-3829

BGH reverse priming site: bases 3874-3891

BGH polyadenylation signal: bases 3880-4104

Kanamycin resistance gene: bases 4277-5071

pUC origin: bases 5371-6044

# Technical Service

---

## World Wide Web



Visit the [Invitrogen Web Resource](#) using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- Subscribe to InvitroWire™, our individually-tailored electronic news service
- View and download vector maps and sequences
- Download manuals in Adobe® Acrobat® (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Post a question at one of our many user forums
- Request catalog and product literature

Once connected to the Internet, launch your web browser (Netscape 3.0 or newer), then enter the following location (or URL):

**<http://www.invitrogen.com>**

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

---

## Phone and E-mail

If you need technical information or help, please e-mail, call, or fax us:

Location	E-mail	Telephone/Fax
United States Canada	tech_service@invitrogen.com	Voice: 1-800-955-6288 Fax: 1-760-602-6500
Mexico South America	tech_service@invitrogen.com	Voice: 01-760-603-7200 Fax: 01-760-602-6500
Europe Africa Middle East	tech_service@invitrogen.nl	Voice: 00800 5345 5345 (free) Fax: 00800 7890 7890 (free) Voice: +31 (0) 50 5299 299 Fax: +31 (0) 50 5299 281
Asia Australia India	pacific_rim@invitrogen.com	Voice: 01-760-603-7200, x250 Fax: 01-760-602-6500

---

## Addresses

If you want to write to us, here are our addresses:

**U.S. Headquarters:**  
Invitrogen Corporation  
1600 Faraday Avenue  
Carlsbad, CA 92008 USA

**European Headquarters:**  
Invitrogen BV  
PO Box 2312, 9704 CH Groningen  
The Netherlands

---

*continued on next page*

## Technical Service, continued

---

### Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, please contact our Technical Service Representatives at 800-955-6288 extension 2 (U.S. and Canada) or 1-760-603-7200 extension 2 (all others).

Invitrogen warrants that all of its products will perform according to the specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.

**Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.**

---

## References

---

- Andersson, S., Davis, D. L., Dahlbäck, H., Jörnvall, H., and Russell, D. W. (1989). Cloning, Structure, and Expression of the Mitochondrial Cytochrome P-450 Sterol 26-Hydroxylase, a Bile Acid Biosynthetic Enzyme. *J. Biol. Chem.* *264*, 8222-8229.
- Ausubel, F. M., Brent, R., Kingston, R. E., Moore, D. D., Seidman, J. G., Smith, J. A., and Struhl, K. (1994). *Current Protocols in Molecular Biology* (New York: Greene Publishing Associates and Wiley-Interscience).
- Bolivar, F., Rodriguez, R. L., Betlach, M. C., and Boyer, H. W. (1977). Construction and Characterization of New Cloning Vehicles. I. Ampicillin-Resistant Derivatives of the Plasmid pMB9. *Gene* *2*, 75-93.
- Bolivar, F., Rodriguez, R. L., Greene, P. J., Betlach, M. C., Heyneker, H. L., and Boyer, H. W. (1977). Construction and Characterization of New Cloning Vehicles. II. A Multipurpose Cloning System. *Gene* *2*, 95-113.
- Boshart, M., Weber, F., Jahn, G., Dorsch-Häsler, K., Fleckenstein, B., and Schaffner, W. (1985). A Very Strong Enhancer is Located Upstream of an Immediate Early Gene of Human Cytomegalovirus. *Cell* *41*, 521-530.
- Chen, C., and Okayama, H. (1987). High-Efficiency Transformation of Mammalian Cells by Plasmid DNA. *Mol. Cell. Biol.* *7*, 2745-2752.
- Chu, G., Hayakawa, H., and Berg, P. (1987). Electroporation for the Efficient Transfection of Mammalian Cells with DNA. *Nuc. Acids Res.* *15*, 1311-1326.
- Davies, J., and Smith, D. I. (1978). Plasmid-Determined Resistance to Antimicrobial Agents. *Ann. Rev. Microbiol.* *32*, 469-518.
- Felgner, P. L., Gadek, T. R., Holm, M., Roman, R., Chan, H. W., Wenz, M., Northrop, J. P., Ringold, G. M., and Danielsen, M. (1987). Lipofectin: A Highly Efficient, Lipid-mediated DNA-transfection Procedure. *Proc. Natl. Acad. Sci. USA* *84*, 7413-7417.
- Felgner, P. L., Holm, M., and Chan, H. (1989). Cationic Liposome Mediated Transfection. *Proc. West. Pharmacol. Soc.* *32*, 115-121.
- Goodwin, E. C., and Rottman, F. M. (1992). The 3'-Flanking Sequence of the Bovine Growth Hormone Gene Contains Novel Elements Required for Efficient and Accurate Polyadenylation. *J. Biol. Chem.* *267*, 16330-16334.
- Kozak, M. (1987). An Analysis of 5'-Noncoding Sequences from 699 Vertebrate Messenger RNAs. *Nuc. Acids Res.* *15*, 8125-8148.
- Kozak, M. (1991). An Analysis of Vertebrate mRNA Sequences: Intimations of Translational Control. *J. Cell Biol.* *115*, 887-903.
- Kozak, M. (1990). Downstream Secondary Structure Facilitates Recognition of Initiator Codons by Eukaryotic Ribosomes. *Proc. Natl. Acad. Sci. USA* *87*, 8301-8305.
- Miller, J. H. (1972). *Experiments in Molecular Genetics* (Cold Spring Harbor, New York: Cold Spring Harbor Laboratory).
- Nelson, J. A., Reynolds-Kohler, C., and Smith, B. A. (1987). Negative and Positive Regulation by a Short Segment in the 5'-Flanking Region of the Human Cytomegalovirus Major Immediate-Early Gene. *Mol. Cell. Biol.* *7*, 4125-4129.

---

*continued on next page*



## References, continued

---

Sambrook, J., Fritsch, E. F., and Maniatis, T. (1989). *Molecular Cloning: A Laboratory Manual, Second Edition* (Plainview, New York: Cold Spring Harbor Laboratory Press).

Shigekawa, K., and Dower, W. J. (1988). Electroporation of Eukaryotes and Prokaryotes: A General Approach to the Introduction of Macromolecules into Cells. *BioTechniques* 6, 742-751.

Wigler, M., Silverstein, S., Lee, L.-S., Pellicer, A., Cheng, Y.-C., and Axel, R. (1977). Transfer of Purified Herpes Virus Thymidine Kinase Gene to Cultured Mouse Cells. *Cell* 11, 223-232.

---

©1998-2001 Invitrogen Corporation. Reproduction forbidden without permission.